

# Headspace Analysis of Volatile Contaminants in Water

Volatile contaminants in waste or waste water from industry can pose a threat to aquatic systems. Examples are the use of chlorinated solvents for cleaning purposes, the accidental release of fuel compounds and the contamination of water by substances released from contaminated industrial sites.

Monitoring of these compounds is required at  $\mu\text{g/L}$  quantities for the protection of human health and aquatic ecosystems. As releases of these compounds are almost impossible to remediate, it is of great importance to investigate rising trends in environmental concentrations well before critical concentrations are reached.

A headspace gas chromatography mass spectrometry system has been set-up in order to determine volatile compounds in water. After equilibration at elevated temperature a sample from the headspace is injected into a GC with Cooled Injection System (CIS). The mass spectrometric detection allows identification of contaminants and the isotope dilution quantification with deuterated or  $\text{C}_{13}$  internal standards. The system is highly automated and reduces the risk of cross-contamination to a minimum. The Headspace-GC/MS consists of an Agilent Technologies 6890 gas chromatograph with 5973N mass spectrometer, two GERSTEL CIS injection systems, a GERSTEL MultiPurpose Sampler (MPS 2) and a GERSTEL Thermal Desorption System (TDS A).

**Table 1:** Volatile compounds which were examined with headspace gas chromatography and detected with mass spectrometry.

Compounds	QS for In-Tr waters (WFD)*	LoD-HS-SIM	LoQ-HS-SIM
Dichloromethane	8.2 $\mu\text{g/L}$	0.010 $\mu\text{g/L}$	0.050 $\mu\text{g/L}$
Trichloromethane	3.85 $\mu\text{g/L}$	0.100 $\mu\text{g/L}$	0.250 $\mu\text{g/L}$
Carbon tetrachloride	7.2 $\mu\text{g/L}$	0.500 $\mu\text{g/L}$	1 $\mu\text{g/L}$
Benzene	16a $\mu\text{g/L}$	0.100 $\mu\text{g/L}$	0.200 $\mu\text{g/L}$
Trichloroethylene	10 $\mu\text{g/L}$	0.002 $\mu\text{g/L}$	0.005 $\mu\text{g/L}$
Tetrachloroethylene	10 $\mu\text{g/L}$	0.001 $\mu\text{g/L}$	0.005 $\mu\text{g/L}$
1,2,4 Trichlorobenzene	1.8 $\mu\text{g/L}$	0.005 $\mu\text{g/L}$	0.010 $\mu\text{g/L}$
1,2,3 Trichlorobenzene		0.005 $\mu\text{g/L}$	0.010 $\mu\text{g/L}$
Naphthalene	2.4 $\mu\text{g/L}$	0.001 $\mu\text{g/L}$	0.005 $\mu\text{g/L}$
Hexachlorobutadiene	0.003 $\mu\text{g/L}$	0.002 $\mu\text{g/L}$	0.005 $\mu\text{g/L}$
1,2 Dichloroethane	1060 $\mu\text{g/L}$	0.500 $\mu\text{g/L}$	1 $\mu\text{g/L}$

\* Quality standard referring to the annual average concentration for inland and transitional waters as proposed by the final outcome of the study: "Towards the Derivation of Quality Standards for Priority Substances in the Context of the Water Framework Directive", Lepper P. (2002), Fraunhofer-Institute Molecular Biology and Applied Ecology


## Authors

A. Cantù, R. Tilio, E. Canuti, G. Hanke & S.J. Eisenreich, Institute for Environment and Sustainability, European Commission - Joint Research Centre, I-21020 Ispra (VA), Italy


Inland and Marine Waters Unit, Organic Analytical Laboratory.

The method meets required performance criteria for drinking water analysis and allows even much lower detection, down to few  $\text{ng/L}$  for analysis of concentration trends. The method can be applied to a minimum concentration of 0.1  $\mu\text{g/L}$  in SCAN-mode (compound dependent) allowing target and non-target analysis while for target analysis in SIM-mode (single ion monitoring) the method was successfully tested with a minimum concentration of 0.005  $\mu\text{g/L}$  (compound dependent).

**Figure 1:** 50  $\text{ng/L}$  hexachlorobutadiene measured in SCAN-mode. Ion 225 is used to quantitate also in sim mode (Figure 2).




**Figure 2:** 50  $\text{ng/L}$  hexachlorobutadiene measured in SIM-mode. It is evident that peaks are better defined and sharper.




# Analysis of 60 volatile compounds with headspace GC / MS

- 1) ?
- 2) Dichlorofluoromethane
- 3) Chloromethane
- 4) Vinyl chloride
- 5) Bromomethane
- 6) Chloroethane
- 7) Trichlorofluoromethane
- 8) Bromodichloromethane
- 9) cis - 1,2 - Dichloroethylene
- 10) t rans - 1,2 - Dichloroethylene
- 11) Methylene chloride
- 12) 1,1 - Dichloroethylene
- 13) tert -Butylmethylether
- 14) 1,1 - Dichloroethane
- 15) 2,2 - Dichloropropane
- 16) Chloroform
- 17) 1,1 - Dichloropropylene
- 18) 1,2 - Dichloroethane
- 19) Trichloroethylene
- 20) Bromochloromethane
- 21) 1,1,1 - Trichloroethane
- 22) Carbon tetrachloride
- 23) 1,2 - Dichloropropane
- 24) Dibromomethane
- 25) Benzene
- 26) Toluene
- 28) Tetrachloroethylene
- 29) Bromoform
- 30) Chlorobenzene
- 31) Styrene
- 32) p-Xylene
- 33) o-Xylene
- 35) 1,2 - Dibromomethane
- 36) 1,3 - Dichloropropane
- 37) 1,1,2 - Trichloroethane
- 38) t rans - 1,3 - Dichloropropylene
- 39) cis - 1,3 - Dichloropropylene
- 40) 1,1,1,2 - Tetrachloroethane
- 41) 1,1,2,2 - Tetrachloroethane
- 42) Ethylbenzene+m-Xylene
- 43) Dibromochloromethane
- 44) 1,2,3 - Trichloropropane
- 46) Isopropylbenzene
- 47) n-Propylbenzene
- 48) 2-Chlorotoluene
- 49) 4-Chlorotoluene
- 50) tert -Butylbenzene
- 51) sec-Butylbenzene
- 52) 1,3-Dichlorobenzene
- 53) 1,4-Dichlorobenzene
- 54) 1,2-Dichlorobenzene
- 55) Bromobenzene
- 56) 1,3,5-Trimethylbenzene
- 57) 1,2,4-Trimethylbenzene
- 58) p-Isopropyltoluene
- 59) n-Butylbenzene
- 60) 1,2,4-Trichlorobenzene
- 61) Naphthalene
- 62) 1,2,3-Trichlorobenzene
- 63) 1,2-Dibromo-3-chloropropane
- 64) Hexachlorobutadiene



Currently 60 compounds are targeted



Identification by comparison with library

## Instrument parameters

### MultiPurpose Sampler (MPS 2) (GERSTEL, autosampler)

Sample volume 1.5 mL  
Incubation temperature 60 °C  
Incubation time 20.00 min  
Syringe temp. 70 °C

### Cooled Injection System (CIS 4 PLUS) (GERSTEL)

Equilibration time 0.2 min  
Initial temperature - 40 °C  
Initial time 0.15 min  
Rate 12 C / min  
Final temperature 200 °C  
Final time 3 min

### 5973 GC / MSD: HS / GC (GC 6890 Series, Agilent Technologies)

Iso time SB: 0 min  
Iso temperature SB: 40 °C  
Iso time 1: 10 min  
Rate: 8 °C / min  
Iso temperature 2: 200 °C  
Iso time 2: 5 min  
Zone temperature transfer line 230 °C  
Capillary column: DB-VRX (J&W)  
Column length: 30 m  
Column i.d.: 0.32 mm  
Film thickness: 1.8 µm  
Glass vaporization tube filled with Tenax TA

### HP 5973 GC / MSD: MSD 5973

Mass range 35 – 500  
Solvent delay 1.00 min  
Zone temperature MS Source 230 °C  
Zone temperature Quad 150 °C

## GERSTEL distribution news

GERSTEL is represented in more than 40 countries world-wide by our own sales and support organization and by a network of distributors. Our distributors are carefully selected according to their experience with analytical instrumentation and their ability to represent and support GERSTEL modules and solutions. Support is a critical factor in offering the best customer value as part of the overall solution. Our distributors can rely on the GERSTEL International Service Support at all times. Only recently, the list of countries with GERSTEL distributors has been extended to include **Egypt, Bahrain, Kuwait, Oman and Qatar**. For contact details, please see the distributor list in this magazine.

In order to provide the best possible support for our many customers in **The Netherlands**, GERSTEL products will be represented exclusively by *Da Vinci*

*Europe BV* as of October 1, 2003. Centrally located in Schiedam, the Netherlands, *Da Vinci Europe BV* offers extensive experience in selling and supporting analytical instruments and in servicing Agilent products.

We are sorry to report the loss of our friend and long-time representative in **Poland**, Mr. Marek Daniewski Ph.D., from the company Chromtec. Our condolences and best wishes go out to Marek's family.

The new GERSTEL distributor in Poland is the company *Perlan Technologies Polska Sp. z o.o.* based in Warsaw. *Perlan Technologies* was formed from the former Agilent Technologies organization in Poland and is an authorized Agilent Technologies distributor and service agent. We are certain that *Perlan Technologies* will provide outstanding support to the Polish GERSTEL customers.

