

# Automated Hydrolysis, Extraction and Analysis of Synthetic Cathinones in Urine using a Robotic Autosampler and LC-MS/MS Platform

Fred D. Foster<sup>1</sup> and John R. Stuff<sup>1</sup>

*1. GERSTEL, Inc., 701 Digital Drive, Suite J, Linthicum, MD 21090, United States*

## KEYWORDS

Sample Preparation, LC-MS/MS, High Throughput Lab Automation, DPX, Bath Salts, Urine, forensic

## ABSTRACT

Synthetic cathinones (also known as Bath Salts) are a group of drug compounds designated as Novel Psychoactive Substances (NPS). They are unregulated, mind-altering substances with no actual approved medicinal use. Since they are cheap substitutes for other stimulants like methamphetamine and cocaine, users will unfortunately turn to these addictive and in some cases even more dangerous alternatives to achieve the desired euphoric effects.

There is a critical need for forensic, health care, and law enforcement scientists to be able to quickly assess and monitor which synthetic cathinone is involved, in order to effectively respond to cases involving these compounds.

Automating the entire hydrolysis, extraction, and subsequent analysis by LC-MS/MS provides the critically needed high throughput analysis for synthetic cathinones in urine. Using the GERSTEL MPS robotic autosampler, syringe transfer of all liquids involved in the enzymatic hydrolysis procedure, controlled mixing of the samples for a defined period, as well as extractions of the subsequent hydrolyzed urine samples using dispersive solid phase extraction were performed. The resulting eluents from the automated extractions were then introduced into the Agilent Ultivo LC-MS/MS instrument.

## INTRODUCTION

A variety of sample handling steps are required prior to the analysis of urine samples to accurately determine analyte concentrations. These steps typically begin with the enzymatic hydrolysis of analytes from their conjugated forms to the native drug using enzymes such as beta-glucuronidase. To ensure that the hydrolysis process is complete and reproducible, the pH and time of the hydrolysis must be controlled and optimized for the enzyme used. The enhanced recombinant beta-glucuronidase B-one from Kura Biotech Inc., is optimized to perform instant hydrolysis at room temperature. The B-One master mix is designed for high-throughput laboratories that use automation and can hydrolyze multiple drug classes at room temperature within 5 minutes with high efficiency<sup>1</sup>.

As a result of this study, we were able to show that an automated enzymatic hydrolysis and subsequent cleanup method was successfully implemented using the GERSTEL MPS robotic sampler for synthetic cathinones in urine. Synthetic cathinone analytes isolated from hydrolyzed urine samples using the automated cleanup procedure were introduced to the LC-MS/MS system, an Agilent Technologies 1260 HPLC coupled with an Agilent Ultivo Triple Quadrupole Mass Spectrometer with Jet stream electrospray source. Combining dispersive solid phase extraction with LC-MS/MS minimizes matrix interference from these biological samples.

**EXPERIMENTAL***Materials:*

All stock solutions for the compounds listed in Table 1 were purchased from Cerilliant. An intermediate analyte stock solution was prepared by combining the analyte stock solution with acetonitrile, resulting in appropriate concentrations for the method evaluation for the different synthetic cathinones.

Compound Name	Precursor Ion (m/Z)	Product Ion (m/z)		Frag (V)		CE (V)		Ret Time (min.)	Ret Time (min.)	Limit of Quantitation (ng/mL)
D3-morphine <sup>1</sup>	289.0	165.1	152.0	153	153	38	66	0.72	-	-
Morphine <sup>1</sup>	286.2	165.1	152.0	158	158	45	64	0.72	-	-
D8-MDPV <sup>2</sup>	284.2	149.0	134.6	120	120	25	25	3.70	-	-
MDPV	276.2	165.0	126.1	120	120	15	20	3.70	1000	2.50
D3-Mephedrone <sup>3</sup>	181.1	163.1	148.0	100	100	5	15	3.23	-	-
Mephedrone <sup>3</sup>	178.1	160.0	144.9	100	100	5	15	3.23	1000	2.50

1 - d<sub>3</sub>-Morphine used as internal standard

2 - d<sub>8</sub>-MDPV used as internal standard

3 - d<sub>3</sub>-Mephedrone used as internal standard

*Table 1. Mass Spectrometer Acquisition Parameters.*

Deuterated analogues, d<sub>3</sub>-mephedrone and d<sub>8</sub>-MDPV were purchased from Cerilliant. An internal standard stock solution containing the deuterated internal standards was prepared in methanol at a concentration of 16.0 µg/mL.

High concentration calibration standard and intermediate QC urine samples were prepared by making appropriate dilutions of the combined intermediate analyte stock solution using analyte-free urine to reach the concentrations listed in Table 1. Calibration standards were then prepared using a dilution ratio strategy from the high concentration sample of 1:2:2:2.5:2:2.5:2:2. The high, middle, and low QC samples were prepared using a dilution ratio strategy from the high concentration sample of 1:10:10. Table 1 lists the concentrations for the highest calibration standard and the limit of quantitation found during analyses.

The enhanced, recombinant β-Glucuronidase **B-One** (cat. #B-One-50mL) was provided from Kura Biotech, Inc. Fresh urine was obtained from a male volunteer. Morphine-6β-D-glucuronide, morphine, and d<sub>3</sub>-morphine were purchased from Cerilliant. All other reagents and solvents used were reagent grade.

*Instrumentation:*

All automated Prep Sequences were performed using a MPS robotic<sup>PRO</sup> sampler with GERSTEL DPX and Agitator Options as shown in Figure 1. All analyses were performed using an Agilent 1260 HPLC with an Agilent Poroshell 120 EC-C18 column, (3.0 x 50mm, 2.7 µm) and an Agilent Ultivo Triple Quadrupole Mass Spectrometer with Jet stream electrospray source. Sample injections were made using the GERSTEL LCMS Tool into a 6 port (0.25mm) Cheminert C2V injection valve outfitted with a 2µL stainless steel sample loop.



*Figure 1. MPS robotic<sup>PRO</sup> sampler with GERSTEL DPX Option.*

### Urine Sample Pretreatment:

1. Pipette 100  $\mu\text{L}$  of urine sample into a clean 2mL autosampler vial and cap with a magnetically transportable cap.

### Automated Prep Sequence:

#### Hydrolysis

1. The MPS adds 20  $\mu\text{L}$  of the internal standard solution to the urine sample.
2. The MPS adds 200  $\mu\text{L}$  of the B-One solution to the urine sample.
3. The vial is moved to the Agitator Option where the sample is mixed at room temperature for 15 minutes while mixing at 750 rpm.
4. The MPS transfers 250  $\mu\text{L}$  of the hydrolyzed urine sample into a clean, empty shell vial.

Figure 2 shows a graphical representation of the general DPX extraction process. The automated DPX extraction used for this method consisted of the following steps:

#### DPX Extraction

1. The MPS aspirates 750 $\mu\text{L}$  of acetonitrile into the 2.5 mL syringe.
2. The MPS picks up a new DPX-RP-S tip.
3. The MPS dispenses 750  $\mu\text{L}$  of acetonitrile through the DPX tip, into the hydrolyzed urine sample located in the shell vial.
4. The MPS aspirates the entire sample followed by 1400  $\mu\text{L}$  of air into the DPX tip to dispersively mix the sample and sorbent.
5. After equilibrating for 5 seconds, the MPS dispenses the contents back into the shell vial and the disposes the DPX tip into a waste container.
6. The MPS transfers 900  $\mu\text{L}$  of 0.1 % formic acid in water into a new, capped 2 mL auto-sampler vial.
7. The MPS finally transfers 100  $\mu\text{L}$  of the supernatant from the DPX extracted hydrolyzed urine sample into the same 2 mL autosampler vial.




Figure 2. Graphical representation of the DPX extraction process.

#### LC Method Parameters

Mobile Phase: A – 0.1% formic acid in water  
B – 0.1% formic acid in acetonitrile

#### LC Gradient:

Time (min)	Flow (mL/min)	Pressure (bar)	% B
0	0.3	800	5
1.0	0.3	800	5
3.0	0.3	800	50
4.0	0.3	800	95
6.0	0.3	800	95
6.1	0.3	800	5

Run time: 10 minutes.  
Injection volume: 2.0 $\mu\text{L}$  (loop over-fill technique)  
Column Temperature: 60  $^{\circ}\text{C}$

#### Mass Spectrometer Parameters:

Operation: Electrospray positive mode  
Gas Temperature: 350  $^{\circ}\text{C}$   
Gas Flow (N<sub>2</sub>): 5 L/min  
Nebulizer pressure: 35 psi  
Sheath Gas Flow (N<sub>2</sub>): 11 L/min  
Sheath Gas Temperature: 400  $^{\circ}\text{C}$   
Capillary voltage: 4000 V  
Nozzle voltage: 500 V  
Delta EMV: 0 V

The mass spectrometer acquisition parameters are shown in Table 1 with qualifier ions.

## RESULTS AND DISCUSSION

To evaluate the automated hydrolysis step of the method, a 150 ng/mL morphine-6β-D-glucuronide sample was prepared in urine. Triplicate 100 μL aliquots of the 150 ng/mL morphine-6β-D-glucuronide urine sample were then hydrolyzed for 0, 0.1, 5, 10, 15, 20, and 30 minutes, respectively, and extracted and analyzed for morphine using the DPX-LC-MS/MS method. As shown in the graph in Figure 3, morphine reached a maximum response after 15 minutes of incubation, proving that morphine had been deconjugated from the glucuronide. An incubation time of 15 minutes was finally chosen to ensure complete and reproducible room temperature hydrolysis of all synthetic cathinones using the B-One solution.




Figure 3. Results from evaluation of room temperature hydrolysis time for morphine using B-One.

Figure 4 shows representative mass chromatograms for the synthetic cathinones, along with their respective qualifier ion transitions, resulting from an extracted low QC sample.




Figure 4. Mass Chromatogram Overlay for hydrolyzed, extracted low QC sample.

The lower limits of quantitation for this method are shown in Table 1. Representative calibration curves are shown in Figures 5 A and B. Regression analysis for both mephedrone and MDPV resulted in R<sup>2</sup> values of 0.99 or greater.




Figure 5 A-B. Representative Calibration Curves: Mephedrone (upper) and MDPV (lower).

The accuracy and precision of the method were determined for the synthetic cathinones analyzed using QC samples at high, middle, and low concentrations. Table 2 shows the resulting accuracy and precision data for the drug compounds. Accuracy data averaged 99.4% (range: 95.7% -102%) and precision data averaged 6.04% CV (range: 2.59% -19.3%) for the synthetic cathinone compounds analyzed.

Compound	QC Level	Exp. Conc. [ng/mL]	Avg. Conc. [ng/mL]	Avg. Precision [%]	Avg. Accuracy [%]
Mephedrone	low	7.50	7.69	5.94	102
	middle	75.0	72.4	3.92	96.5
	high	750	743	3.48	99.1
MDPV	low	7.50	7.17	19.3	95.7
	middle	75.0	74.4	3.87	99.2
	high	750	777	2.59	104

Table 2. QC sample accuracy and precision table.

## CONCLUSIONS

As a result of this study, we were able to show:

- Synthetic cathinones in urine samples can be successfully hydrolyzed, extracted using an automated SPE procedure, and determined using the Agilent Ultivo Triple Quadrupole Mass Spectrometer.
- This method was readily automated using the GERSTEL MPS robotic<sup>PRO</sup> sampler.
- Evaluation of the hydrolysis method showed morphine to be deconjugated within 15 minutes at room temperature.
- Linear calibration curves resulting in R<sup>2</sup> values 0.99 or greater were achieved for the determined synthetic cathinones.
- The hydrolysis-SPE-LC-MS/MS method proved to be accurate and precise. Accuracy data averaged 99.4% (range: 95.7% -102%) and precision data averaged 6.04 %CV (range: 2.59% -19.3%) for the determined synthetic cathinone compounds.

## REFERENCES

- (1) “Technical Data Sheet, Product: B-One”, Kura Biotech Inc., Retrieved January 2020 from <http://www.kurabiotec.com/wp-content/uploads/2019/12/Technical-DataSheet-B-One.pdf>

**GERSTEL GmbH & Co. KG**

Eberhard-Gerstel-Platz 1  
45473 Mülheim an der Ruhr  
Germany

+49 (0) 208 - 7 65 03-0  
+49 (0) 208 - 7 65 03 33  
gerstel@gerstel.com  
www.gerstel.com

---

## GERSTEL Worldwide

**GERSTEL, Inc.**

701 Digital Drive, Suite J  
Linthicum, MD 21090  
USA

+1 (410) 247 5885  
+1 (410) 247 5887  
sales@gerstelus.com  
www.gerstelus.com

**GERSTEL AG**

Wassergrabe 27  
CH-6210 Sursee  
Switzerland

+41 (41) 9 21 97 23  
+41 (41) 9 21 97 25  
swiss@ch.gerstel.com  
www.gerstel.ch

**GERSTEL K.K.**

1-3-1 Nakane, Meguro-ku  
Tokyo 152-0031  
SMBC Toritsudai Ekimae Bldg 4F  
Japan

+81 3 5731 5321  
+81 3 5731 5322  
info@gerstel.co.jp  
www.gerstel.co.jp

**GERSTEL LLP**

10 Science Park Road  
#02-18 The Alpha  
Singapore 117684

+65 6779 0933  
+65 6779 0938  
SEA@gerstel.com  
www.gerstel.com

**GERSTEL (Shanghai) Co. Ltd**

Room 206, 2F, Bldg.56  
No.1000, Jinhai Road,  
Pudong District  
Shanghai 201206

+86 21 50 93 30 57  
china@gerstel.com  
www.gerstel.cn

**GERSTEL Brasil**

Av. Pascoal da Rocha Falcão, 367  
04785-000 São Paulo - SP Brasil

+55 (11)5665-8931  
+55 (11)5666-9084  
gerstel-brasil@gerstel.com  
www.gerstel.com.br

Information, descriptions and specifications in this Publication are subject to change without notice. GERSTEL, GRAPHPACK and TWISTER are registered trademarks of GERSTEL GmbH & Co. KG.

© Copyright by GERSTEL GmbH & Co. KG



Awarded for the active pursuit of environmental sustainability